

## CYP4F2 Conjugated Antibody

Catalog No: #C44314



Package Size: #C44314 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CYP4F2 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	CYP4F2 Antibody detects endogenous levels of total CYP4F2
Immunogen Description	The antiserum was produced against synthesized peptide derived from human CYP4F2.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CYP4F2, CYP4F2, Cytochrome P450 4F2, Cytochrome P450-LTB-omega, CPF2
Accession No.	Swiss-Prot#:P78329NCBI Gene ID:8529NCBI mRNA#:NCBI Protein#:
Calculated MW	60
Formulation	Sodium Phosphate, 0.5% BSA, 40% glycerin, 0.02% sodium azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.