

## ZNF497 Antibody

Catalog No: #44746



Package Size: #44746-1 50ul #44746-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF497 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu
Specificity	ZNF497antibody detects endogenous levels of total?ZNF497
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZNF497
Target Name	ZNF497
Other Names	Zinc finger prothein 497
Accession No.	Swiss-Prot#: Q6ZNH5NCBI Gene ID:162968
Calculated MW	55 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western blot analysis of ZNF497 expression in HEK293 cells ,The lane on the left is treated with the antigen-specific peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.