

JAM1 Antibody

Catalog No: #44829



Package Size: #44829-1 50ul #44829-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JAM1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	JAM1 antibody detects endogenous levels of total JAM1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human JAM1
Target Name	JAM1
Other Names	CD 321 antibody CD321 antibody CD321 antibody CD321 antigen antibody ESTM33 antibody F11 receptor antibody F11R antibody JAM 1 antibody JAM A antibody JAM antibody JAM-1 antibody JAM-A antibody JAM1 antibody JAM1_HUMAN antibody JAMA antibody JCAM antibody Jcam1 antibody Junction adhesion molecule 1
Accession No.	Swiss-Prot#: Q9Y624NCBI Gene ID:50848
Calculated MW	33 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.