

## CPD Conjugated Antibody

Catalog No: #C45850



Package Size: #C45850 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CPD Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Human
Specificity	CPD Antibody detects endogenous levels of total CPD
Immunogen Description	A synthesized peptide derived from human CPD
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CPD, Carboxypeptidase D, gp180, Metallo-carboxypeptidase D
Accession No.	Swiss-Prot#:O75976NCBI Gene ID:1362NCBI mRNA#:NCBI Protein#:
Calculated MW	153
Formulation	Sodium Phosphate, 0.5% BSA, 40% glycerin, 0.02% sodium azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.