

LUC7L3 Conjugated Antibody

Catalog No: #C45886

Package Size: #C45886 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LUC7L3 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	HumanO PredictedO Rat,Ms
Specificity	LUC7L3 Antibody detects endogenous levels of total LUC7L3
Immunogen Description	Recombinant fragment of human LUC7L3 protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	LUC7L3, CRA, CRE-associated protein 1, CREAP-1, CROP, hLuc7A, Luc7-like protein 3, LUC7-like 3 (S. cerevisiae), LUC7A, O48, OA48-18, CREAP1
Accession No.	Swiss-Prot#:O95232NCBI Gene ID:51747NCBI mRNA#:NCBI Protein#:
Calculated MW	51
Formulation	Sodium Phosphate, 0.5% BSA, 40% glycerin, 0.02% sodium azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.