

LGR5 Antibody

Catalog No: #44974



Package Size: #44974-1 50ul #44974-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	LGR5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB, IH, IF/IC
Species Reactivity	Hu Ms Rt
Specificity	LGR5antibody detects endogenous levels of total?LGR5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LGR5
Target Name	LGR5
Other Names	FEX antibody G protein coupled receptor 49 antibody G protein coupled receptor 67 antibody g protein-coupled receptor fex antibody G-protein coupled receptor 49 antibody G-protein coupled receptor 67 antibody G-protein coupled receptor HG38 antibody GPR 49 antibody
Accession No.	Swiss-Prot#: O75473NCBI Gene ID:8549
Calculated MW	100 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.