

JAK3 (Phospho-Tyr785) Polyclonal Conjugated Antibody

Catalog No: #C12331

Orders: order@signalwayantibody.com
 Package Size: #C12331-AF350 100ul #C12331-AF405 100ul #C12331-AF488 100ul #C12331-AF555 100ul #C12331-AF594 100ul

 Support: tech@signalwayantibody.com

#C12331-AF647 100ul #C12331-AF680 100ul #C12331-AF750 100ul #C12331-Biotin 100ul #C12331-Conjugated 50ul

Description

Product Name	JAK3 (Phospho-Tyr785) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu
Specificity	Phospho-JAK3 (Y785) Polyclonal Antibody detects endogenous levels of JAK3 protein only when phosphorylated at Y785.
Immunogen Description	Synthesized peptide derived from human JAK3 around the phosphorylation site of Y785.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	JAK3; Tyrosine-protein kinase JAK3; Janus kinase 3; JAK-3; Leukocyte janus kinase; L-JAK
Accession No.	Swiss-Prot#:P52333NCBI Gene ID:3718
Calculated MW	125
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.