

JNK1/2/3 (Phospho-Tyr185) Polyclonal Conjugated Antibody

Catalog No: #C12345

Orders: order@signalwayantibody.com
 Package Size: #C12345-AF350 100ul #C12345-AF405 100ul #C12345-AF488 100ul #C12345-AF555 100ul #C12345-AF594 100ul

 Support: tech@signalwayantibody.com

#C12345-AF647 100ul #C12345-AF680 100ul #C12345-AF750 100ul #C12345-Biotin 100ul #C12345-Conjugated 50ul

Description

Product Name	JNK1/2/3 (Phospho-Tyr185) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu
Specificity	Phospho-JNK1/2/3 (Y185) Polyclonal Antibody detects endogenous levels of JNK1/2/3 protein only when phosphorylated at Y185.
Immunogen Description	Synthesized peptide derived from human JNK1/2/3 around the phosphorylation site of Y185.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MAPK8; JNK1; PRKM8; SAPK1; SAPK1C; Mitogen-activated protein kinase 8; MAP kinase 8; MAPK 8; JNK-46; Stress-activated protein kinase 1c; SAPK1c; Stress-activated protein kinase JNK1; c-Jun N-terminal kinase 1; MAPK9; JNK2; PRKM9; SAPK1A; Mi
Accession No.	Swiss-Prot#:P45983/P45984/P53779NCBI Gene ID:5599/5601/5602
Calculated MW	45
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.