

## JARID1C Conjugated Monoclonal Antibody

Catalog No: #C27002



Package Size: #C27002-AF350 100ul #C27002-AF405 100ul #C27002-AF488 100ul #C27002-AF555 100ul #C27002-AF594 100ul  
 #C27002-AF647 100ul #C27002-AF680 100ul #C27002-AF750 100ul #C27002-Biotin 100ul #C27002-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JARID1C Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Applications	WB, IF
Species Reactivity	Hu Ms Rt D
Immunogen Description	Recombinant fragment(Human)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	lysine-specific demethylase 5C, protein SmcX, selected cDNA on X, Smcy homolog, X-linked, Smcx homolog, X chromosome, histone demethylase JARID1C, JmjC domain-containing protein SMCX, Jumonji/ARID domain-containing protein 1C
Accession No.	Swiss-Prot#:P41229 NCBI Gene ID:8242
Calculated MW	220
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.