

VDAC1 antibody

Catalog No: #23076

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	VDAC1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 181 of VDAC1
Target Name	VDAC1
Accession No.	NCBI Gene ID: 7416NCBI mRNA#: BC008482NCBI Protein#: AAH08482
Concentration	0.4mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 31kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:500

Images




Sample (30 ug of whole cell lysate)

A: A431

12% SDS PAGE

VDAC1 antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Cal27 xenograft, using VDAC1 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.