

MPG Conjugated Antibody

Catalog No: #C32837

Package Size: #C32837-AF350 100ul #C32837-AF405 100ul #C32837-AF488 100ul #C32837-AF555 100ul #C32837-AF594 100ul #C32837-AF647 100ul #C32837-AF680 100ul #C32837-AF750 100ul #C32837-Biotin 100ul #C32837-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MPG Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB,IHC,ELISA,IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total MPG protein.
Immunogen Description	Recombinant protein of human MPG.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AAG;MDG;ADPG;APNG;Mid1
Accession No.	Swiss-Prot#:P29372NCBI Gene ID:4350
Calculated MW	32
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Antibodies were purified by affinity purification using immunogen.

Background

Hydrolysis of the deoxyribose N-glycosidic bond to excise 3-methyladenine, and 7-methylguanine from the damaged DNA polymer formed by alkylation lesions.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.