

FAP-1 Conjugated Antibody

Catalog No: #C33369

Package Size: #C33369-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FAP-1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Peptide affinity purification
Applications	WB, IF
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total FAP-1 protein.
Immunogen Description	Synthesized peptide derived from human FAP-1.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DPPIV;FAPA;SEPRASE;fibroblast activation protein alpha subunit;integral membrane serine protease
Accession No.	Swiss-Prot#:Q12884NCBI Gene ID:2191
Calculated MW	88
Formulation	Phosphate buffered saline with 0.5% BSA, 0.02% sodium azide and 40% Glycerol
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Background

In association with DPP4 is involved in the pericellular proteolysis of the extracellular matrix (ECM), the migration and invasion of endothelial cells into the ECM. May have a role in tissue remodeling during development and wound healing, and may contribute to invasiveness in malignant cancers.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.