

JAK1 Conjugated Antibody

Catalog No: #C33721

Package Size: #C33721-AF350 100ul #C33721-AF405 100ul #C33721-AF488 100ul #C33721-AF555 100ul #C33721-AF594 100ul #C33721-AF647 100ul #C33721-AF680 100ul #C33721-AF750 100ul #C33721-Biotin 100ul #C33721-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JAK1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IF
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total JAK1 protein.
Immunogen Description	Synthesized peptide derived from internal of human JAK1.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	EC 2.7.10.2;JAK-1;JAK1;JAK1A;Janus kinase 1
Accession No.	Swiss-Prot#:P23458NCBI Gene ID:3716
Calculated MW	130
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

Product Description

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Background

Tyrosine kinase of the non-receptor type, involved in the IFN-alpha/beta/gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.