

TAF9B Antibody

Catalog No: #45238



Package Size: #45238-1 50ul #45238-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAF9B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms
Specificity	TAF9B Antibody detects endogenous levels of total TAF9B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TAF9B
Target Name	TAF9B
Other Names	TAF9B, DN7, TAFII31L, TAF9L, TFIID-31, DN-7, TBP-associated factor 9L
Accession No.	Swiss-Prot#: Q9HBM6NCBI Gene ID:51616
Calculated MW	28 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis TAF9B using NIH-3T3 whole cell lysates

Published Papers

el at., miR-146a attenuates apoptosis and modulates autophagy by targeting TAF9b/P53 pathway In doxorubicin-Induced cardiotoxicity.

In Cell Death Dis on 2019 Sep 11 by Pan J, Tang Y et al..PMID:31511497, , (2019)

[PMID:31511497](https://pubmed.ncbi.nlm.nih.gov/31511497/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.