

## ZBTB17 Antibody

Catalog No: #45519



Package Size: #45519-1 50ul #45519-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZBTB17 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ZBTB17 Antibody detects endogenous levels of total ZBTB17
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZBTB17
Target Name	ZBTB17
Other Names	ZBTB17, MIZ-1, MIZ1, Zinc finger protein 151, ZNF60, PHZ-67, Zinc finger protein 60, ZNF151
Accession No.	Swiss-Prot#: Q13105NCBI Gene ID:7709
Calculated MW	88 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis ZBTB17 using HeLa whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.