

JDP2 Polyclonal Conjugated Antibody

Catalog No: #C41083

Package Size: #C41083-AF350 100ul #C41083-AF405 100ul #C41083-AF488 100ul #C41083-AF555 100ul #C41083-AF594 100ul #C41083-AF647 100ul #C41083-AF680 100ul #C41083-AF750 100ul #C41083-Biotin 100ul #C41083-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JDP2 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu Ms Rt
Specificity	JDP2 Polyclonal Antibody detects endogenous levels of JDP2 protein.
Immunogen Description	Synthesized peptide derived from human JDP2 around the non-phosphorylation site of T148.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	JDP2; Jun dimerization protein 2
Accession No.	Swiss-Prot#:Q8WYK2NCBI Gene ID:122953
Calculated MW	22
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.