

JIP-3 Polyclonal Conjugated Antibody

Catalog No: #C41084



Package Size: #C41084-AF350 100ul #C41084-AF405 100ul #C41084-AF488 100ul #C41084-AF555 100ul #C41084-AF594 100ul #C41084-AF647 100ul #C41084-AF680 100ul #C41084-AF750 100ul #C41084-Biotin 100ul #C41084-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JIP-3 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu Ms
Specificity	JIP-3 Polyclonal Antibody detects endogenous levels of JIP-3 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human JIP-3.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MAPK8IP3; JIP3; KIAA1066; C-Jun-amino-terminal kinase-interacting protein 3; JIP-3; JNK-interacting protein 3; JNK MAP kinase scaffold protein 3; Mitogen-activated protein kinase 8-interacting protein 3
Accession No.	Swiss-Prot#:Q9UPT6NCBI Gene ID:23162
Calculated MW	147
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.