

JNK3 Polyclonal Conjugated Antibody

Catalog No: #C41087

Package Size: #C41087-AF350 100ul #C41087-AF405 100ul #C41087-AF488 100ul #C41087-AF555 100ul #C41087-AF594 100ul #C41087-AF647 100ul #C41087-AF680 100ul #C41087-AF750 100ul #C41087-Biotin 100ul #C41087-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JNK3 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu Ms Rt
Specificity	JNK3 Polyclonal Antibody detects endogenous levels of JNK3 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human JNK3.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MAPK10; JNK3; JNK3A; PRKM10; SAPK1B; Mitogen-activated protein kinase 10; MAP kinase 10; MAPK 10; MAP kinase p49 3F12; Stress-activated protein kinase 1b; SAPK1b; Stress-activated protein kinase JNK3; c-Jun N-terminal kinase 3
Accession No.	Swiss-Prot#:P53779NCBI Gene ID:5602
Calculated MW	48
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.