

Oct-3/4 Polyclonal Conjugated Antibody

Catalog No: #C41262



Package Size: #C41262-AF350 100ul #C41262-AF405 100ul #C41262-AF488 100ul #C41262-AF555 100ul #C41262-AF594 100ul #C41262-AF647 100ul #C41262-AF680 100ul #C41262-AF750 100ul #C41262-Biotin 100ul #C41262-Conjugated 50ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Oct-3/4 Polyclonal Conjugated Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Applications | WB, IF |
| Species Reactivity | Hu Ms Rt |
| Specificity | Oct-3/4 Polyclonal Antibody detects endogenous levels of Oct-3/4 protein. |
| Immunogen Description | Synthesized peptide derived from the Internal region of human Oct-3/4. |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | POU5F1; OCT3; OCT4; OTF3; POU domain; class 5, transcription factor 1; Octamer-binding protein 3; Oct-3; Octamer-binding protein 4; Oct-4; Octamer-binding transcription factor 3; OTF-3 |
| Accession No. | Swiss-Prot#:Q01860NCBI Gene ID:5460 |
| Calculated MW | 50 |
| Formulation | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide |
| Storage | Store at 4°C in dark for 6 months |

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.