

Op18 Polyclonal Conjugated Antibody

Catalog No: #C41287

Package Size: #C41287-AF350 100ul #C41287-AF405 100ul #C41287-AF488 100ul #C41287-AF555 100ul #C41287-AF594 100ul
 #C41287-AF647 100ul #C41287-AF680 100ul #C41287-AF750 100ul #C41287-Biotin 100ul #C41287-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Op18 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu Ms Rt
Specificity	Op18 Polyclonal Antibody detects endogenous levels of Op18 protein.
Immunogen Description	Synthesized peptide derived from human Op18 around the non-phosphorylation site of S25.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	STMN1; C1orf215; LAP18; OP18; Stathmin; Leukemia-associated phosphoprotein p18; Metablastin; Oncoprotein 18; Op18; Phosphoprotein p19; pp19; Prosolin; Protein Pr22; pp17
Accession No.	Swiss-Prot#:P16949NCBI Gene ID:3925
Calculated MW	19
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.