

JAM3 Antibody

Catalog No: #45639



Package Size: #45639-1 50ul #45639-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JAM3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	JAM3 Antibody detects endogenous levels of total JAM3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human JAM3
Target Name	JAM3
Other Names	JAM3, JAMC, Junctional adhesion molecule C, JAM-2, Junctional adhesion molecule 3, JAM-3, JAM-C
Accession No.	Swiss-Prot#: Q9BX67NCBI Gene ID:83700
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.