

Transferrin Polyclonal Conjugated Antibody

Catalog No: #C41645



Package Size: #C41645-Conjugated 50ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transferrin Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB;IHC;IF;ELISA
Species Reactivity	Hu Ms Rt
Specificity	Transferrin Polyclonal Antibody detects endogenous levels of Transferrin protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Transferrin.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TF; Serotransferrin; Transferrin; Beta-1 metal-binding globulin; Siderophilin
Accession No.	Swiss-Prot#:P02787NCBI Gene ID:7018
Calculated MW	77
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB 1:500 - 1:2000.

IHC-p: 1:100-300

ELISA: 1:20000..

IF 1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.