

SLAM Polyclonal Conjugated Antibody

Catalog No: #C41685



Package Size: #C41685-AF350 100ul #C41685-AF405 100ul #C41685-AF488 100ul #C41685-AF555 100ul #C41685-AF594 100ul #C41685-AF647 100ul #C41685-AF680 100ul #C41685-AF750 100ul #C41685-Biotin 100ul #C41685-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	SLAM Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu Rt Ms
Specificity	SLAM Polyclonal Antibody detects endogenous levels of SLAM protein.
Immunogen Description	Synthesized peptide derived from the region of human SLAM.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	SLAMF1; SLAM; Signaling lymphocytic activation molecule; CDw150; IPO-3; CD150
Accession No.	Swiss-Prot#:Q13291NCBI Gene ID:6504
Calculated MW	37
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.