

Jagged1 Polyclonal Conjugated Antibody

Catalog No: #C41786

Package Size: #C41786-AF350 100ul #C41786-AF405 100ul #C41786-AF488 100ul #C41786-AF555 100ul #C41786-AF594 100ul #C41786-AF647 100ul #C41786-AF680 100ul #C41786-AF750 100ul #C41786-Biotin 100ul #C41786-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Jagged1 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu Ms Rt
Specificity	Jagged1 Polyclonal Antibody detects endogenous levels of Jagged1 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Jagged1.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	JAG1; JAGL1; Protein jagged-1; Jagged1; hJ1; CD339
Accession No.	Swiss-Prot#:P78504NCBI Gene ID:182
Calculated MW	140
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.