

T3E Polyclonal Conjugated Antibody

Catalog No: #C41832



Package Size: #C41832-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T3E Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WBoO IFoO IHCoo IF
Species Reactivity	HuoO RatoO MouseoO Monkey
Specificity	T3E Polyclonal Antibody detects endogenous levels of T3E protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human T3E.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD3E; T3E; T-cell surface glycoprotein CD3 epsilon chain; T-cell surface antigen T3/Leu-4 epsilon chain; CD3e;CD3- α
Accession No.	Swiss-Prot#:P07766NCBI Gene ID:916
Calculated MW	23
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.