

JAM-A Polyclonal Conjugated Antibody

Catalog No: #C41850

Package Size: #C41850-AF350 100ul #C41850-AF405 100ul #C41850-AF488 100ul #C41850-AF555 100ul #C41850-AF594 100ul #C41850-AF647 100ul #C41850-AF680 100ul #C41850-AF750 100ul #C41850-Biotin 100ul #C41850-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JAM-A Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Hu Rt
Specificity	JAM-A Polyclonal Antibody detects endogenous levels of JAM-A protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human JAM-A.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	F11R; JAM1; JCAM; Junctional adhesion molecule A; JAM-A; Junctional adhesion molecule 1; JAM-1; Platelet F11 receptor; Platelet adhesion molecule 1; PAM-1; CD321
Accession No.	Swiss-Prot#:Q9Y624NCBI Gene ID:50848
Calculated MW	32
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.