

## SOX17 Antibody

Catalog No: #45703



Package Size: #45703-1 50ul #45703-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	SOX17 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB,IHC
Species Reactivity	Hu Ms
Specificity	SOX17 Antibody detects endogenous levels of total SOX17
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SOX17
Target Name	SOX17
Other Names	SOX17, Transcription factor SOX-17, VUR3
Accession No.	Swiss-Prot#: Q9H6I2NCBI Gene ID:64321
Calculated MW	47 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details


WB 1:500-2000;

IHC 1:50-300

## Images



Western Blot analysis of extracts from Jurkat cells, using Sox-17 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human liver cancer.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.