

# Vasopressin V1b receptor Polyclonal ;Conjugated Antibody

Catalog No: #C42598

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #C42598-AF350 100ul #C42598-AF405 100ul #C42598-AF488 100ul #C42598-AF555 100ul #C42598-AF594 100ul

#C42598-AF647 100ul #C42598-AF680 100ul #C42598-AF750 100ul #C42598-Biotin 100ul #C42598-Conjugated 50ul

## Description

Product Name	Vasopressin V1b receptor Polyclonal ;Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IF
Species Reactivity	Rt
Specificity	The antibody detects endogenous level of total Vasopressin V1b receptor polyclonal antibody;.
Immunogen Description	Recombinant Rat Vasopressin V1b receptor protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AVPR V1b,AVPR V3,Vasopressin V3 receptor
Accession No.	Swiss-Prot#:P48974
Calculated MW	47
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Background

Receptor for arginine vasopressin. The activity of this receptor is mediated by G proteins which activate a phosphatidyl-inositol-calcium second messenger system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.