

YY1AP1 Conjugated Antibody

Catalog No: #C42866

Package Size: #C42866-AF350 100ul #C42866-AF405 100ul #C42866-AF488 100ul #C42866-AF555 100ul #C42866-AF594 100ul #C42866-AF647 100ul #C42866-AF680 100ul #C42866-AF750 100ul #C42866-Biotin 100ul #C42866-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	YY1AP1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total YY1AP1 protein.
Immunogen Description	Fusion protein of human YY1AP1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	YAP; HCCA1; HCCA2; YY1AP
Accession No.	Swiss-Prot#:Q9H869NCBI Gene ID:55249NCBI mRNA#:BC009358
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

Background

The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms.?

Note: This product is for in vitro research use only and is not intended for use in humans or animals.