

YPEL1 Conjugated Antibody

Catalog No: #C43586

Package Size: #C43586-AF350 100ul #C43586-AF405 100ul #C43586-AF488 100ul #C43586-AF555 100ul #C43586-AF594 100ul #C43586-AF647 100ul #C43586-AF680 100ul #C43586-AF750 100ul #C43586-Biotin 100ul #C43586-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	YPEL1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total YPEL1 protein.
Immunogen Description	Full length fusion protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	FKSG3
Accession No.	Swiss-Prot#:O60688NCBI Gene ID:29799NCBI Protein#:BC074501
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

Background

This gene is located in the region associated with DiGeorge syndrome on chromosome 22. The encoded protein localizes to the centrosome and nucleolus and may play a role in the regulation of cell division.?

Note: This product is for in vitro research use only and is not intended for use in humans or animals.