

CLM7 Polyclonal Conjugated Antibody

Catalog No: #C46830

Package Size: #C46830-AF350 100ul #C46830-AF405 100ul #C46830-AF488 100ul #C46830-AF555 100ul #C46830-AF594 100ul #C46830-AF647 100ul #C46830-AF680 100ul #C46830-AF750 100ul #C46830-Biotin 100ul #C46830-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CLM7 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IHC;IF;ELISA
Species Reactivity	Hu;Rat;Mouse
Specificity	The antibody detects endogenous CLM7
Immunogen Description	Synthetic peptide from human CLM7 protein at AA range: 61-110
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CMRF35-like molecule 7 (CLM-7) (CD300 antigen-like family member B) (CMRF35-A2) (Immune receptor expressed on myeloid cells 3) (IREM-3) (Leukocyte mono-Ig-like receptor 5) (Triggering receptor expressed on myeloid cells 5) (TREM-5; CD300b)
Accession No.	Swiss-Prot#:A8K4G0NCBI Gene ID:124599
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Immunohistochemistry: IHC 1:50-300

ELISA 1:10000-20000.

IF 1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.