

JPT2 Conjugated Antibody

Catalog No: #C47286

Package Size: #C47286-AF350 100ul #C47286-AF405 100ul #C47286-AF488 100ul #C47286-AF555 100ul #C47286-AF594 100ul #C47286-AF647 100ul #C47286-AF680 100ul #C47286-AF750 100ul #C47286-Biotin 100ul #C47286-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JPT2 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total JPT2 protein.
Immunogen Description	Fusion protein of human JPT2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	L11; HN1L; C16orf34
Accession No.	Swiss-Prot#:Q9H910NCBI Gene ID:90861NCBI Protein#:BC014438
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

IF: 1:50-1:200

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.