

JAZF1 Conjugated Antibody

Catalog No: #C47673

Package Size: #C47673-AF350 100ul #C47673-AF405 100ul #C47673-AF488 100ul #C47673-AF555 100ul #C47673-AF594 100ul #C47673-AF647 100ul #C47673-AF680 100ul #C47673-AF750 100ul #C47673-Biotin 100ul #C47673-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JAZF1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IF
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total JAZF1 protein.
Immunogen Description	Synthetic peptide of human JAZF1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TIP27; ZNF802
Accession No.	Swiss-Prot#:Q86VZ6NCBI Gene ID:221895NCBI Protein#:NP_778231
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

IF: 1:50-1:200

Background

This gene encodes a nuclear protein with three C2H2-type zinc fingers, and functions as a transcriptional repressor. Chromosomal aberrations involving this gene are associated with endometrial stromal tumors. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.