

JMY Conjugated Antibody

Catalog No: #C47894



Package Size: #C47894-AF350 100ul #C47894-AF405 100ul #C47894-AF488 100ul #C47894-AF555 100ul #C47894-AF594 100ul #C47894-AF647 100ul #C47894-AF680 100ul #C47894-AF750 100ul #C47894-Biotin 100ul #C47894-Conjugated 50ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	JMY Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IHC
Species Reactivity	Hu,Ms
Specificity	Recognizes endogenous levels of JMY protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human JMY.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Junction-mediating and -regulatory protein
Accession No.	Swiss-Prot#:Q8N9B5NCBI Gene ID:133746
Calculated MW	133
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.