

VPS72 Conjugated Antibody

Catalog No: #C47943

Package Size: #C47943-AF350 100ul #C47943-AF405 100ul #C47943-AF488 100ul #C47943-AF555 100ul #C47943-AF594 100ul
 #C47943-AF647 100ul #C47943-AF680 100ul #C47943-AF750 100ul #C47943-Biotin 100ul #C47943-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VPS72 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB, IHC, IF/ICC
Species Reactivity	Hu, Ms
Specificity	Recognizes endogenous levels of VPS72 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human VPS72.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TCFL1; YL1; Vacuolar protein sorting-associated protein 72 homolog; Protein YL-1; Transcription factor-like 1
Accession No.	Swiss-Prot#:Q15906NCBI Gene ID:6944
Calculated MW	41
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.