

Vitronectin Conjugated Antibody

Catalog No: #C48876



Package Size: #C48876-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vitronectin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	12D1
Applications	WB, IF
Species Reactivity	Hu, Ms, Rt
Immunogen Description	recombinant protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Complement S Protein antibody Epibolin antibody S Protein antibody S-protein antibody Serum Spreading Factor antibody Serum-spreading factor antibody Somatomedin B antibody Somatomedin-B antibody V75 antibody Vitronectin V10 subunit antibody Vitronectin V65 subunit antibody VN antibody VNT antibody VTN antibody VTNC_HUMAN antibody
Accession No.	Swiss-Prot#:P04004
Calculated MW	75 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Background

A431, A549, MCF-7, HepG2, mouse liver tissue, mouse skin tissue, human liver tissue, human liver cancer tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.