

LRP4 Antibody

Catalog No: #46034

Package Size: #46034-1 50ul #46034-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRP4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	LRP4 Antibody detects endogenous levels of total LRP4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LRP4
Target Name	LRP4
Other Names	LRP4, CLSS, LRP-4, SOST2, KIAA0816, MEGF7
Accession No.	Swiss-Prot#: O75096NCBI Gene ID:4038
Calculated MW	212kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis LRP4 using NIH-3T3 whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.