

DCPS Polyclonal Antibody

Catalog No: #27360

Package Size: #27360-1 50ul #27360-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DCPS Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human DCPS (NP_054745.1).
Other Names	DCPS; ARS; DCS1; HINT-5; HINT5; HSL1; HSPC015; m7GpppX diphosphatase
Accession No.	Swiss-Prot#:Q96C86NCBI Gene ID:28960
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:1000 - 1:2000

Images



Western blot analysis of extracts of HT-29 cells, using DCPS antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.