

## DNAJB2 Polyclonal Antibody

Catalog No: #27389

Package Size: #27389-1 50ul #27389-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DNAJB2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human DNAJB2 (NP_006727.2).
Other Names	DNAJB2; CMT2T; DSMA5; HSJ-1; HSJ1; HSPF3; dnaJ homolog subfamily B member 2
Accession No.	Swiss-Prot#:P25686NCBI Gene ID:3300
Calculated MW	42kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:200 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using DNAJB2 antibody.

## Background

This gene is almost exclusively expressed in the brain, mainly in the neuronal layers. It encodes a protein that shows sequence similarity to bacterial DnaJ protein and the yeast homologs. In bacteria, this protein is implicated in protein folding and protein complex dissociation. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.