

ZC3H14 Polyclonal Antibody

Catalog No: #27399

Package Size: #27399-1 50ul #27399-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZC3H14 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZC3H14 (NP_997545.2).
Other Names	ZC3H14; MRT56; MSUT-2; NY-REN-37; SUT2; UKp68; zinc finger CCCH-type containing 14
Accession No.	Swiss-Prot#:Q6PJT7NCBI Gene ID:79882
Calculated MW	105kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using ZC3H14 antibody.

Background

The protein encoded by this gene is a poly(A)-binding protein that can affect gene expression and poly(A) tail length. The encoded protein may influence mRNA stability, nuclear export, and translation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.