

ENTPD3 Polyclonal Antibody

Catalog No: #27430

Package Size: #27430-1 50ul #27430-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ENTPD3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human ENTPD3 (NP_001278890.1).
Other Names	ENTPD3; CD39L3; HB6; NTPDase-3; ectonucleoside triphosphate diphosphohydrolase 3
Accession No.	Swiss-Prot#:O75355NCBI Gene ID:956
Calculated MW	59kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse lung, using ENTPD3 antibody at 1:1000 dilution.

Background

This gene encodes a plasma membrane-bound divalent cation-dependent E-type nucleotidase. The encoded protein is involved in the regulation of extracellular levels of ATP by hydrolysis of it and other nucleotides. Multiple transcript variants have been described.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.