

# JOSD1 Polyclonal Antibody

Catalog No: #27432



Package Size: #27432-1 50ul #27432-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	JOSD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human JOSD1 (NP_055691.1).
Other Names	JOSD1; dJ508115.2; josephin-1
Accession No.	Swiss-Prot#:Q15040NCBI Gene ID:9929
Calculated MW	27kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using JOSD1 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.