

ZDHHC18 Antibody

Catalog No: #46120



Package Size: #46120-1 50ul #46120-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZDHHC18 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ZDHHC18 Antibody detects endogenous levels of total ZDHHC18
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZDHHC18
Target Name	ZDHHC18
Other Names	ZDHHC18, DHHC-18, DHHC18, Palmitoyltransferase ZDHHC18
Accession No.	Swiss-Prot#: Q9NUE0NCBI Gene ID:84243
Calculated MW	42 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis ZDHHC18 using LOVO whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.