

ZSCAN4 Polyclonal Antibody

Catalog No: #27590

Package Size: #27590-1 50ul #27590-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZSCAN4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZSCAN4 (NP_689890.1).
Other Names	ZSCAN4; ZNF494; zinc finger and SCAN domain containing 4
Accession No.	Swiss-Prot#:Q8NAM6NCBI Gene ID:201516
Calculated MW	49kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of rat brain, using ZSCAN4 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.