

DESI2 Polyclonal Antibody

Catalog No: #27623

Package Size: #27623-1 50ul #27623-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DESI2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human DESI2 (NP_057160.2).
Other Names	DESI2; C1orf121; CGI-146; DESI; DESI1; DeSI-2; FAM152A; PNAS-4; PPPDE1; desumoylating isopeptidase 2
Accession No.	Swiss-Prot#:Q9BSY9NCBI Gene ID:51029
Calculated MW	21kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DESI2 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.