

## TUBB8 Polyclonal Antibody

Catalog No: #27733

Package Size: #27733-1 50ul #27733-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TUBB8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human TUBB8
Other Names	TUBB8;OOMD;OOMD2;bA631M21.2
Accession No.	Uniprot:Q3ZCM7GeneID:347688
Calculated MW	54kDa
SDS-PAGE MW	55kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


## Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using TUBB8 antibody.



Immunofluorescence analysis of U2OS cells using TUBB8 antibody.

## Background

The protein encoded by this gene represents the primary beta-tubulin subunit of oocytes and the early embryo. Defects in this gene, which is primate-specific, are a cause of oocyte maturation defect 2 and infertility.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.