

EPB41L3 Polyclonal Antibody

Catalog No: #27774

Package Size: #27774-1 50ul #27774-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EPB41L3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human EPB41L3 (NP_001268463.1).
Other Names	EPB41L3; 4.1B; DAL-1; DAL1; band 4.1-like protein 3
Accession No.	Swiss-Prot#:Q9Y2J2NCBI Gene ID:23136
Calculated MW	121kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using EPB41L3 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.