

ZBTB21 Polyclonal Antibody

Catalog No: #27784

Package Size: #27784-1 50ul #27784-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZBTB21 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZBTB21 (NP_065778.3).
Other Names	ZBTB21; ZNF295; zinc finger and BTB domain containing 21
Accession No.	Swiss-Prot#:Q9ULJ3NCBI Gene ID:49854
Calculated MW	120kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZBTB21 antibody at 1:3000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.