

# NKAIN4 Polyclonal Antibody

Catalog No: #27810



Package Size: #27810-1 50ul #27810-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NKAIN4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NKAIN4 (NP_690603.3).
Other Names	NKAIN4; C20orf58; FAM77A; bA261N11.2; sodium/potassium transporting ATPase interacting 4
Accession No.	Swiss-Prot#:Q8IVV8NCBI Gene ID:128414
Calculated MW	Refer to Figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:1000 - 1:2000

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.