

## DLGAP1 Polyclonal Antibody

Catalog No: #27877



Package Size: #27877-1 50ul #27877-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DLGAP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human DLGAP1 (NP_001229690.1).
Other Names	DLGAP1; DAP-1; DAP-1-ALPHA; DAP-1-BETA; DAP1; DLGAP1A; DLGAP1B; GKAP; SAPAP1; DLG associated protein 1
Accession No.	Swiss-Prot#:O14490NCBI Gene ID:9229
Calculated MW	109kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using DLGAP1 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.